1

**Instructions:**

Warning: As with any exercise program, consult a health care professional before using the Handmaster Plus system. If pain is experienced, discontinue use immediately and consult a health care professional. Do not overuse.

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* Tip: To load the cord into the ball, use an unwrapped paperclip (not included). Select the large or small cord for hand size and resistance desired.

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Use placement tabs to avoid twisting and damage to elastic loops. The loops fit tight by design.

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
**Finger and Thumb Placement:**

Place the thumb loop on the base of the thumb first. Place the finger loops between the top two knuckles on each of the fingers. Place the loops on the correct finger. (thumb=1, small finger=5, etc.). The Handmaster Plus can be used on the right hand (RH) or the left hand (LH).

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**Finger Thumb Hand and Carpal Tunnel:**

Squeeze the ball for one second. The ball provides resistance. Open and spread the fingers and thumb for one second. The cord provides resistance. Repeat until comfortable fatigue only.

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**Wrist Forearm and Elbow :**

Squeeze the ball for one second. Extend only the fingers for one second. Extend the wrist for one second. Repeat until comfortable fatigue only.

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
**Fingers and Thumb:**

Squeeze the ball, using the fingertips only, for one second. Open and spread the fingers and thumb for one second. Repeat until comfortable fatigue only.

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